

# T4 Polynucleotide Kinase

**Catalogue number:** MB00801, 500 U

## Description

T4 Polynucleotide Kinase (T4 PNK) catalyses the transfer of the terminal phosphate of ATP to 5'-hydroxyl termini of polynucleotides such as DNA and RNA, oligonucleotides and 3'-mononucleotides. This enzyme also possesses a 3'-phosphatase and 2', 3' cyclic phosphodiesterase activity. The enzyme can be used for the following applications: 5'-phosphorylation of oligonucleotides, PCR products, other DNA or RNA to allow subsequent ligation; end-labelling DNA or RNA for use as probes, for DNA sequencing or for DNA-protein foot printing; removal of 3' phosphoryl groups.

## Storage temperature

T4 Polynucleotide Kinase should be stored at -20 °C in a constant temperature freezer. The protein will remain stable till the expiry date if stored as specified.

## Unit definition

One enzyme unit incorporates 1 nmol of radiolabelled ATP into DNA substrate in 30 min at 37 °C under standard assay conditions.

**Enzyme concentration:** 30 U/μL

## Inactivation

T4 Polynucleotide Kinase is heat inactivated at 65 °C for 10 min.

## System components and Reaction conditions

T4 Polynucleotide Kinase is provided with a dedicated highly optimized NZYTech reaction buffer, as well as a dilution buffer to optimize levels of incorporation. See section "Important notes" below for further information. The optimum temperature for activity is 37 °C.

## End labelling protocol (radioactive 5'-end labelling)

This protocol is suitable for 10 pmol of dephosphorylated DNA or oligonucleotide.

1. Prepare the following 50 μL reaction:

Oligonucleotide (10 pmol)	x μL
[ $\gamma$ - <sup>32</sup> P]ATP*	y μL
10× Reaction buffer for T4 PNK	5 μL
Nuclease-free H <sub>2</sub> O	Up to 50 μL
T4 Polynucleotide Kinase (diluted)	2-5 units
Final volume	50 μL

\* The number of pmol of [ $\gamma$ -<sup>32</sup>P]ATP should be at least 2 times the number of pmol of 5'-ends.

2. Mix gently and pulse.

3. Incubate at 37 °C for 30 min.

4. Stop the reaction by heating at 65 °C for 10 minutes. Proceed with the separation of 5'-end-labeled oligonucleotide from precursor ATP by thin-layer or column chromatography according to standard methods.

## DNA phosphorylation protocol

The following standard protocol serves as a general guideline for a non-radioactive kinase reaction. Preferably the enzyme should be added last.

1. Prepare the following 50 μL reaction:

Linear dsDNA or oligonucleotide*	x μL
ATP (10 mM) (not provided)	5 μL
10× Reaction buffer for T4 PNK	5 μL
Nuclease-free H <sub>2</sub> O	Up to 50 μL
T4 Polynucleotide Kinase (diluted)	10 units
Final volume	50 μL

\* up to 300 pmol of 5-end.

2. Mix gently and pulse.

3. Incubate at 37 °C for 30 min.

4. Stop the reaction by heating at 65 °C for 10 minutes.

## Important notes

- Dilute the T4 Polynucleotide Kinase using the Dilution buffer provided. Adding concentrated enzyme can result in lower levels of incorporation. Do not store the enzyme in the diluted form.
- Labelling efficiency is influenced by the terminal 5' nucleotide. An oligonucleotide with G at the 5'-end labels about 6-fold higher than a 5'C and about 1.5-fold higher than a 5'T or 5'A.

## Quality control assays

### Purity

T4 Polynucleotide Kinase is >95% pure as judged by SDS polyacrylamide gel electrophoresis followed by BlueSafe (NZYTech, Cat. No. MB152) staining.

### Nuclease assays

To test for DNase contamination, 0.2-0.3 μg of supercoiled pNZY28 plasmid DNA are incubated with 10 U of T4 Polynucleotide Kinase for 14-16 hours at 37 °C. To test for RNase contamination, 1 μg of RNA is incubated with 10 U of T4 Polynucleotide Kinase for 1 hour at 37 °C. Following incubation, the nucleic acids are visualized on a GreenSafe-stained agarose gel. There must be no visible nicking or cutting of the nucleic acids.

### Functional assay

T4 polynucleotide Kinase is tested for activity in a reaction to phosphorylate an insert DNA required for ligation with a non-phosphorylated vector. The efficiency of Kinase's reaction is evaluated by the number of bacterial colonies transformed with the cloning product.

## Related products

Product name	Cat. No.
Water for Molecular Biology	MB11101
Alkaline Phosphatase	MB018
Speedy Ligase	MB130
NZYGelpure	MB011

## Certificate of Analysis

Test	Result
Enzyme purity	Pass
Nucleases contamination	Pass
Functional assay	Pass

Approved by:



Patricia Ponte  
Senior Manager, Quality Systems

*For research use only.*

