

## T5 Exonuclease

**Catalogue number:** MB43301, 1000 U

### Description

T5 Exonuclease is a typical DNA nuclease that functions in the 5' to 3' direction and attacks both single and double-stranded DNA. The enzyme substrate is the 5'-termini of a DNA molecule in the context of single or double-stranded nucleic acids as well as gaps and nicks of linear or circular double-stranded DNA. However, T5 Exonuclease does not degrade supercoiled double-stranded DNA. Thus, the enzyme may be used to remove denatured DNA generated during the alkaline-lysis protocols for plasmid purification or linear and circular DNA from supercoiled plasmid preps. It may also be used to remove incomplete ligation products from ligated circular double-stranded DNA.

### Storage conditions

T5 Exonuclease should be stored at -20 °C in a constant temperature freezer. The protein will remain stable till the expiry date if stored as specified.

### Unit definition

One unit of enzyme activity is defined as the amount of enzyme required to cause the change of 0.00032 A260 nm/min at 37 °C in the enzyme dedicated buffer in a total reaction volume of 50 µL.

**Enzyme concentration:** 10 U/µL

### Inactivation

T5 Exonuclease is highly resistant to heat inactivation. Thus, alternative protocols should be considered when requiring removing the enzyme from reactions, such as DNA silica column purification or phenol/chloroform extraction.

### System components and Reaction conditions

T5 Exonuclease is provided with a dedicated and highly optimized NZYtech reaction buffer and displays an optimum temperature of 37 °C.

### Standard protocol

The following standard protocol serves as a general guideline to degrade single-stranded or double-stranded nicked or linear DNA with T5 Exonuclease (the enzyme does not attack supercoiled dsDNA). Preferably the enzyme should be added last.

1. Prepare the following 50 µL reaction:

Component	Volume
Substrate DNA	≤ 1 µg
T5 Exonuclease reaction buffer (10x)	5 µL
T5 Exonuclease	1 µL (10 U)
Nuclease-free H <sub>2</sub> O (Cat. No. MB11101)	up to 50 µL

**Note:** It may be required to titrate the enzyme or test different incubation periods for more specific results or partial digestions.

- Gently mix and pulse.
- Incubate at 37 °C for 30 minutes.
- If required, stop the reaction by adding EDTA to at least 15 mM final concentration.
- To obtain a highly pure product, perform a column purification step using NZYGelpure kit (NZYtech, Cat. No. MB011). Best results may be achieved by separating cleaved DNA through agarose gel electrophoresis prior to DNA clean-up.

### Quality Control Assays

#### Purity

T5 Exonuclease is >95% pure as judged by SDS polyacrylamide gel electrophoresis followed by BlueSafe staining (NZYtech, Cat. No. MB15201).

#### Nucleases assays

To test for DNase contamination, 0.2-0.3 µg of supercoiled pNZY28 plasmid DNA are incubated with 10 U of T5 Exonuclease for 14-16 hours at 37 °C. To test for RNase contamination, 1 µg of RNA is incubated with 10 U of T5 Exonuclease for 1 hour at 37 °C. Following incubation, the nucleic acids are visualized on a GreenSafe-stained agarose gel. There must be no visible nicking or cutting of the nucleic acids.

#### Functional assay

T5 Exonuclease is assayed in a reaction containing 100 ng of a double stranded DNA fragment of 1 kb under standard conditions. Activity is measured by monitoring reduction of A260 over time.

V2401

*For research use only.*