

Mycoplasma Cleaner

Catalogue number	Presentation
MB46401	250 mL
MB46402	500 mL (Refill)

Features

- Ready-to-use solution specifically designed to decontaminate laboratory surfaces and apparatus.
- It offers a wide spectrum of activity, effectively targeting *Mycoplasma* spp., *S. aureus*, *E. faecalis*, *E. coli*, *P. mirabilis*, *P. aeruginosa*, *C. albicans*, as well as HBV/HIV, Papovavirus, Adenovirus, Poliovirus, and other contaminants.
- The cleaner is conveniently packaged in a user-friendly, lightweight bottle equipped with a trigger spray for easy application.
- Using this cleaner simplifies the decontamination process, saving both time and effort.
- A refill option is available for added convenience.

Description

Mycoplasma Cleaner is a highly effective decontaminating solution with broad-spectrum activity, designed to target not only *Mycoplasma* but also viruses, bacteria, and fungi. This solution provides a safe and reliable approach to decontamination in laboratory settings, particularly on benches, safety cabinets, incubators, cell storage boxes, and liquid nitrogen containers, among other surfaces and equipment.

Mycoplasma species, which are common contaminants in cell cultures, can be challenging to detect using standard microscopy techniques due to their small size. As a result, *Mycoplasma* infections often go unnoticed, leading to cell infection and a range of detrimental effects. Contamination can occur through media spills from contaminated cell cultures during loading and unloading procedures in biosafety cabinets, resulting in surface contamination where *Mycoplasma* can survive for extended periods. To prevent cross-contamination, it is crucial to promptly isolate and eliminate infected cultures, as even low levels of *Mycoplasma* infection can disrupt cell growth, alter metabolism, and protein synthesis. Regular and thorough cleaning and disinfection of all surfaces and equipment are vital measures to combat and prevent *Mycoplasma* contamination in the laboratory.

By utilizing the Mycoplasma Cleaner and implementing proper hygiene practices, you can effectively mitigate the risk of mycoplasma contamination, ensuring the integrity of your cell cultures and the reliability of your research.

As an additional measure for efficient contamination control, we highly recommend conducting routine testing of cell culture materials using our ***Mycoplasma* spp. qPCR Kit (Cat. No. MD05511)**. This kit is specifically designed for the in vitro detection of *Mycoplasma* spp. genomes, providing a fast, simple, and effective method to identify mycoplasma-contaminated cell cultures. The real-time qPCR-based method utilized by the *Mycoplasma* spp. qPCR Kit can be conveniently performed in any standard molecular biology laboratory, offering a valuable tool to ensure the integrity and reliability of your cell culture experiments.

The Mycoplasma Cleaner has undergone extensive testing and validation to ensure optimal performance in diverse laboratory settings. Its effectiveness is scientifically proven and supported by robust evidence, providing you with the assurance you need for critical decontamination requirements.

Storage temperature

Store at 15 °C to 25 °C.

Protocol of usage

Use gloves as personal protective equipment when manipulating this product.

Decontamination of surfaces:

1. Spray Mycoplasma Cleaner directly onto the contaminated surfaces until they are completely moistened.
2. Allow the product to act for 5 minutes.
3. Gently wipe the surface with a paper towel until dry. The Cleaner does not leave any residue.

Note: Prior to use, it is recommended to perform a surface compatibility test as Mycoplasma Cleaner may cause discoloration on sensitive surfaces.

Decontamination of apparatus:

1. Apply Mycoplasma Cleaner onto a paper towel.
2. Vigorously rub all surfaces of the apparatus, ensuring thorough coverage.
3. Thoroughly dry the surface with a clean paper towel.

Please note that while our recommendations provide a robust decontamination solution, we encourage users to verify the suitability of Mycoplasma Cleaner for their specific purposes and applications. Do not add the cleaner to the reaction.

For life science research only. Not for use in diagnostic procedures.

NZYtech Lda. Estrada do Paço do Lumiar, Campus do Lumiar - Edifício E, R/C, 1649-038 Lisboa, Portugal Tel.: +351.213643514 Fax:
+351.217151168 www.NZYtech.com