


pHTP16 Expression Vector

pHTP16 was designed for the cloning and expression of high-levels of recombinant proteins in *Escherichia coli*. Recombinant proteins are expressed in fusion with a *Ruminococcus flavefaciens* cellulosomal protein (CpA), which is able to promote solubility and folding of fusion partners. This vector, included in the portfolio of NZYtech pHTP expression vectors, is part of the NZYEasy Cloning & Expression System. pHTP16 contains two poly-histidine (6xHis) sequences (N- and C-terminal) which allow subsequent recombinant protein purification by immobilized metal ion affinity chromatography (IMAC).

1. Vector Map



pHTP16 Cloning/Expression Region

<i>Nco I</i>	CpA	<i>Nco I</i>	His-Tag
<u>CCATGG</u> GTTGGCGAGTGTACGGCTATATCGAC.804bp..GATGTCAAGGATTACTCAATCGTTG		<u>CCATGG</u> GCAGCAGCCATCATCATCATCACAGCAGCGGC	
MetGlyGlyGluCysHisGlyThrIleAsp.268aa..AspValLysAspTyrSerIleValAlaMetGlySerSer		HisHisHisHisHisHisSerSerGly	
CCTCAGCAAGGCTGAGG / ⌘ / CCTCAGCTTCCGCTGAGGTCCGTCGACAAGCTTGC		<u>CTCGAGCACCACCACCACCACCAC</u> TGAGATCCGGCTGCT	
ProGlnGlnGlyLeuArg / ⌘ / ProGlnLeuProLeuArgSerValAspLysLeuAlaAlaAlaLeuGlu		HisHisHisHisHisHis*	STOP

⌘ Represents the site where the gene will be inserted.

Note 1: For correct expression, inserted gene needs to be in frame with pHTP16 5' gene sequence. Inserts correctly cloned into pHTP16 will maintain reading frames starting on the ATG codon.

Note 2: CpA is a recombinant cellulosomal protein (Cp) that is highly expressed in *E. coli*. This protein has carbohydrate-binding module activity, displaying affinity for β -glycans (xyloglucan, glucomannan, galactomannan and barley β -glucan).

2. Vector Sequence (6145 bp)

