

NZY 1X Tris-Glycine-SDS Buffer Powder

Catalogue number	Presentation
MB46801	1 pouch (for 1 L)
MB46802	5 pouches (for 5 L)

Description

NZY 1X Tris-Glycine-SDS Buffer Powder is used for electrophoresis on NZY Tris-Glycine Precast Gel (MB46601) or Laemmli Tris-HCl gels in Tris-Glycine buffer system. It is convenient and universal for electrophoresis in Tris-Glycine buffer system.

Key Features

- **Reliable:** rigorous quality control to ensure consistent separation of proteins in electrophoresis.
- **Convenient:** premeasured pouches make 1 liter of 1X buffer solution without needing pH adjustment.
- **Fast:** dissolves within minutes and it is ready to use.
- **Stable:** powder packaging suitable for extended storage.

Application

Running buffer for Laemmli Tris-HCl gel electrophoresis.

Contents

0.025 M Tris, 0.192 M glycine, 0.10% SDS

Shipping Conditions

This product must be shipped at room temperature.

Storage Conditions

Store NZY 1X Tris-Glycine-SDS Buffer Powder at 15 °C to 25 °C.

For life science research only. Not for use in diagnostic procedures.