

PDZRN3_1 PDZ Domain, *Homo sapiens*

Catalogue number	Presentation
PD01001	0.25 mg
PD01002	3 x 0.25 mg

Description

PDZRN3_1 PDZ Domain from *Homo sapiens* is a recombinant protein purified from *Escherichia coli*. The protein is provided in 35 mM NaHepes buffer, pH 7.5, 750 mM NaCl, 200 mM imidazol, 3.5 mM CaCl₂ and 25% (v/v) glycerol, at a 0.5 mg/mL concentration. Bulk quantities of this product can be made available upon request. To place an order, simply visit our website. We offer fast and secure shipping worldwide.

Electrophoretic Purity

PDZRN3_1 PDZ Domain purity was evaluated using sodium dodecyl sulfate-polyacrylamide gel electrophoresis (SDS-PAGE), followed by BlueSafe staining (MB15201) (Figure 1).




Figure 1. SDS-PAGE analysis of PDZRN3_1 PDZ Domain was conducted in (Lane 1), employing a 14% polyacrylamide gel. The enzyme exhibits a band corresponding to a molecular weight of approximately 13.12 kDa. Lane M contains a Protein Marker for reference.

Storage temperature

The protein should be stored at -30°C to -15°C in a constant temperature freezer. The protein will remain stable till the expiry date if stored as specified.

Protein Sequence

TKSLTLVLHRDSGSLGFNIIGRPSVDNHDGSSSEGIFVSKIVDSGPAAKEGGLQIHDRIIEVNGRDLRATHDQAVEAFKTAKEPIVVQVLRRT

Number of PDZ in native protein

1

Temperature and pH optima

The protein exhibits optimal activity within a pH of 7 and at a temperature of 36.5 °C.

PDB code

1UHP

Reference

Unpublished

Customer Support

Our dedicated customer support team is always ready to assist you with any questions or technical issues you may have. Reach us via email at info@nzytech.com.

Quality control assay

Protein purity is determined to be ≥90%, as assessed by SDS-PAGE and subsequent BlueSafe staining (MB15201).

For life science research only. Not for use in diagnostic procedures.

NZYtech Lda. Estrada do Paço do Lumiar, Campus do Lumiar - Edifício E, R/C, 1649-038 Lisboa, Portugal Tel.: +351.213643514 Fax: +351.217151168 www.nzytech.com