

SARS-CoV-2 Diagnostic Kit Selection Guide

NZYTech has a diverse catalogue of SARS-CoV-2 RT-PCR kits with CE-IVD certification. All kits are designed to work perfectly with various widely used real-time PCR instruments and allow you to have complete control of your determinations - all real-time PCR kits include a positive control that validates experimental performance and an internal control to confirm efficient nucleic acids extraction from human biological samples, as well as absence of PCR inhibitors.



TARGET PATHOGEN

			Determinations	SARS-CoV-2	Other resp. pathogens	Target Genes	Target Channels	Internal Control (IC)	IC Channel
CE-IVD	MD03201	SARS-CoV-2 One-Step RT-PCR Kit, IVD	96	✔		RdRp	FAM™	Rnase P	HEX™/JOE™/VIC™
CE-IVD	MD04831	SARS-CoV-2 One-Step RT-PCR Kit, RdRp and N target genes, IVD	96	✔		RdRp and N	FAM™	Rnase P	HEX™/JOE™/VIC™
CE-IVD	MD04832	SARS-CoV-2 One-Step RT-PCR Kit, RdRp and N target genes, IVD	4 x 96	✔		RdRp and N	FAM™	Rnase P	HEX™/JOE™/VIC™
CE-IVD	MD04871	SARS-CoV-2 One-Step RT-PCR Kit II, RdRp and N target genes, IVD	96	✔		RdRp and N	HEX™/JOE™/VIC™ and FAM™	Rnase P	Cy5™
CE-IVD	MD04872	SARS-CoV-2 One-Step RT-PCR Kit II, RdRp and N target genes, IVD	4 x 96	✔		RdRp and N	HEX™/JOE™/VIC™ and FAM™	Rnase P	Cy5™
CE-IVD	MD04851	COVID-19 & Flu A/B Multiplex One-Step RT-qPCR Kit, IVD	96	✔	Influenza A Influenza B	RdRp and N / InfA M1 and InfB NS2	TexasRed™ and FAM™	Rnase P	HEX™/JOE™/VIC™
	MD03191	SARS-CoV-2 One-Step RT-PCR Kit	96	✔		RdRp	FAM™	Rnase P	HEX™/JOE™/VIC™
	MD04841	Veterinary SARS-CoV-2 One-Step RT-qPCR Kit	96	✔		N	FAM™	β-actin	HEX™/JOE™/VIC™