

## Laccase A, *Escherichia coli*

### EcLaccA (AA NC)


Catalogue number	Presentation
CZ01621	1 mg
CZ01622	3 x 1 mg

#### Description

Laccase A (EcLaccA), assigned the E.C. number 1.10.3.-, is a derivative of *Escherichia coli*. It is a multicopper oxidase. The recombinant EcLaccA, purified from *Escherichia coli*, is a single-domain Auxiliary Activity enzyme family 0 (AA NC) enzyme (see more details at [www.cazy.org](http://www.cazy.org)). The protein is supplied in a solution containing 35 mM NaHepes buffer (pH 7.5), 750 mM NaCl, 200 mM Imidazole, 3.5 mM CaCl<sub>2</sub>, and 25% (v/v) glycerol, at a concentration of 1 mg/mL. Bulk quantities of this product can be made available upon request. To place an order, simply visit our website. We offer fast and secure shipping worldwide.

#### Electrophoretic Purity

The molecular integrity and purity of EcLaccA (AA NC) were evaluated using sodium dodecyl sulfate-polyacrylamide gel electrophoresis (SDS-PAGE), followed by BlueSafe staining (MB15201) (Figure 1).



**Figure 1.** SDS-PAGE analysis of EcLaccA (AA NC) was conducted in (Lane 1), employing a 14% polyacrylamide gel. The enzyme exhibits a band corresponding to a molecular weight of approximately 55,47 kDa. Lane M contains a Protein Marker for reference.

#### Storage temperature

The protein should be stored at -30°C to -15°C in a constant temperature freezer. The protein will remain stable till the expiry date if stored as specified.

#### Substrate specificity

EcLaccA (AA NC) hydrolyses 2,20-azinobis(3-ethylbenzthiazoline-6-sulfonic acid) (ABTS).

#### Temperature and pH optima

The pH optimum for enzymatic activity is 6.5 while temperature optimum is 55 °C.

## Enzyme activity

The substrate specificity and kinetic properties of *EclaccA* (AA NC) are detailed in the referenced publication provided below. To perform enzyme assays and determine specific activity values, adhere to the methodology outlined in the cited paper(s).

## Reference

Grass and Rensing (2001) *Biochem Biophys Res Commun.* 7;286(5):902-8.

Roberts *et al.* (2002) *Proc Natl Acad Sci U S A.* 99(5):2766-71.

## Customer Support

Our dedicated customer support team is always ready to assist you with any questions or technical issues you may have. Reach us via email at [info@nzytech.com](mailto:info@nzytech.com).

## Quality control assay

Protein purity is determined to be  $\geq 90\%$ , as assessed by SDS-PAGE and subsequent BlueSafe staining (MB15201).

For life science research only. Not for use in diagnostic procedures.