

Cellulase 131A, *Podospora anserina* *PaCel131A* (GH131-CBM1)

Catalogue number:

CZ02961, 0,5 mg
CZ02962, 3 × 0,5 mg

Description

Cellulase 131A (*PaCel131A*), assigned the E.C. number 3.2.1.-, is a derivative of *Podospora anserina*. It is an enzyme that displays a broad exo-1,3/1,6- β -glucanase specificity with endo- β -1,4-glucanase activity. The recombinant *PaCel131A*, purified from *Pichia pastoris*, is a modular Glycoside Hydrolase family 131 (GH131) enzyme, with a C-terminal family 1 Carbohydrate Binding Module (CBM1) (see more details at www.cazy.org). The protein is supplied in a solution containing 35 mM NaHepes buffer (pH 7.5), 750 mM NaCl, 200 mM Imidazole, 3.5 mM CaCl₂, and 25% (v/v) glycerol, at a concentration of 0,5 mg/mL. Bulk quantities of this product can be made available upon request. To place an order, simply visit our website. We offer fast and secure shipping worldwide.

Electrophoretic Purity

The molecular integrity and purity of *PaCel131A* (GH131-CBM1) were evaluated using sodium dodecyl sulfate-polyacrylamide gel electrophoresis (SDS-PAGE), followed by BlueSafe staining (MB15201) (Figure 1).



Figure 1. SDS-PAGE analysis of *PaCel131A* (GH131-CBM1) was conducted in Lane 1, employing a 14% polyacrylamide gel. The enzyme exhibits a singular band corresponding to a molecular weight of approximately 66 kDa. Lane M contains a Protein Marker for reference.

Storage temperature

The protein should be stored at -30°C to -15°C in a constant temperature freezer. The protein will remain stable till the expiry date if stored as specified.

Substrate specificity

PaCel131A (GH131-CBM1) hydrolyses 1,3- β , 1,4- β , 1,3-1,4- β and 1,6- β -glucans displaying highest specificity for laminarin.

Temperature and pH optima

The enzyme exhibits optimal activity within a pH range of 4.5-7.0 and at a temperature of 40°C. Maximal enzymatic activity is achieved at pH 5.5 and a consistent temperature of 40°C.

Enzyme activity

The substrate specificity and kinetic properties of *PaCel131A* (GH131-CBM1) are detailed in the referenced publication provided below. To perform enzyme assays and determine specific activity values, adhere to the methodology outlined in the cited paper(s).

Reference

Lafond *et al.* (2012) Applied and Environmental Microbiology 78, 8540–8546.

Customer Support

Our dedicated customer support team is always ready to assist you with any questions or technical issues you may have. Reach us via email at info@nzytech.com.

Quality control assay

Protein purity is determined to be $\geq 90\%$, as assessed by SDS-PAGE and subsequent BlueSafe staining (MB15201).

V2301

For research use only