

Cellulase 8C, *Clostridium cellulolyticum*

CcCel8C (GH8)

Catalogue number	Presentation
CZ03911	1 mg
CZ03912	3 x 1 mg

Description

Cellulase 8C (CcCel8C), assigned the E.C. number 3.2.1.4, is a derivative of *Clostridium cellulolyticum*. It is an endo-1,4- β -glucanase. The recombinant CcCel8C, purified from *Escherichia coli*, is a single-domain Glycoside Hydrolase family 8 (GH8) enzyme (see more details at www.cazy.org). The protein is supplied in a solution containing 35 mM NaHepes buffer (pH 7.5), 750 mM NaCl, 200 mM Imidazole, 3.5 mM CaCl₂, and 25% (v/v) glycerol, at a concentration of 1 mg/mL. Bulk quantities of this product can be made available upon request. To place an order, simply visit our website. We offer fast and secure shipping worldwide.

Electrophoretic Purity

The molecular integrity and purity of CcCel8C (GH8) were evaluated using sodium dodecyl sulfate-polyacrylamide gel electrophoresis (SDS-PAGE), followed by BlueSafe staining (MB15201) (Figure 1).




Figure 1. SDS-PAGE analysis of CcCel8C (GH8) was conducted in (Lane 1), employing a 14% polyacrylamide gel. The enzyme exhibits a band corresponding to a molecular weight of approximately 42,95 kDa. Lane M contains a Protein Marker for reference.

Storage temperature

The protein should be stored at -30°C to -15°C in a constant temperature freezer. The protein will remain stable till the expiry date if stored as specified.

Substrate specificity

CcCel8C (GH8) hydrolyses 1,4- β -glucans in particular soluble cellulose.

Temperature and pH optima

The pH optimum for enzymatic activity is 6 while temperature optimum is 48 °C.

Enzyme activity

The substrate specificity and kinetic properties of CcCel8C (GH8) are detailed in the referenced publication provided below. To perform enzyme assays and determine specific activity values, adhere to the methodology outlined in the cited paper(s).

Reference

Fierobe *et al.* (1993) *Eur J Biochem.* 217(2):557-65.

Customer Support

Our dedicated customer support team is always ready to assist you with any questions or technical issues you may have. Reach us via email at info@nzytech.com.

Quality control assay

Protein purity is determined to be $\geq 90\%$, as assessed by SDS-PAGE and subsequent BlueSafe staining (MB15201).

For life science research only. Not for use in diagnostic procedures.

NZYtech Lda. Estrada do Paço do Lumiar, Campus do Lumiar - Edifício E, R/C, 1649-038 Lisboa, Portugal Tel.:+351.213643514 Fax:
+351.217151168 www.nzytech.com